

Six-Month and One-Year Clinical Outcomes after Placement of a Dedicated Bifurcation Device:

A Patient-Level Pooled Analysis of Eight Registry Studies Confirms the Predictable Clinical Outcomes¹

Maik J. Grundeken, MD, Solomon Asgedom, MD, Peter Damman, MD, Maciej Lesiak, MD, PhD, Michael S. Norell, MD, Eulogio García, MD, Armando Bethencourt, MD, Pier Woudstra, MD, Karel T. Koch, MD, PhD, Marije M. Vis, MD, Jose P. Henriques, MD, PhD, Yoshinobu Onuma, MD, David P. Foley, MD, PhD, Antonio L. Bartorelli, MD, PhD, Pieter R. Stella, MD, PhD, Jan G. Tijssen, PhD, Robbert J. de Winter, MD, PhD, Joanna J. Wykrzykowska, MD, PhD

Objective of the Analysis:

- Confirm the clinical results from small size registries and studies in a larger sample size
- Evaluate clinical outcomes beyond 6 months

Evaluated Outcomes:

- Primary:
 - 6-month Target Vessel Failure (TVF), composite of cardiac death, any myocardial infarction (MI), clinically indicated target vessel revascularization (TVR)
- Secondary:
 - 1-year TVF
 - 6-month and 1-year cardiac death, any MI, clinically indicated TVR, clinically indicated TLR, Stent Thrombosis

Patient level data merged from **8 registries and studies** included **905 patients** with **929 bifurcations**.

Registry Cohorts	Initiated by	Number of Patients included in the Current Analysis	Follow-up
First-In-Man*	Sponsor	30	6 months
eTryton 150/Benelux**	Sponsor	216 [‡]	6 months
eTryton Spain	Sponsor	142	6 months
IUVANT	Investigator	67	6 months
Amsterdam ⁺	Investigator	91	6 months [†]
Poznan ^{††}	Investigator	91	9 months [†]
Wolverhampton	Investigator	79	18 months [†]
Dublin	Investigator	189	16 months [†]

1. Data presented at EuroPCR 2012

* Onuma Y, Muller R, Ramcharitar S, et al. Tryton I, First-In-Man (FIM) study: six month clinical and angiographic outcome, analysis with new quantitative coronary angiography dedicated for bifurcation lesions. EuroIntervention 2008.

**Agostoni P, Foley D, Lesiak M, et al. A prospective multicentre registry, evaluating real-world usage of the Tryton side branch stent: results of the E-Tryton 150/Benelux registry. EuroIntervention 2012. [‡] In the original publication 302 patients were included; duplicates are removed and reported as part of the Poznan and Dublin cohorts.

[†] Median follow-up.

^{††} Magro M, Wykrzykowska J, Serruys PW, et al. Six-month clinical follow-up of the Tryton side branch stent for the treatment of bifurcation lesions: A two center registry analysis. Catheter Cardiovasc Interv 2011.

⁺ M. Grundeken et al. Six-month clinical outcomes of the Tryton Side Branch Stent for the treatment of bifurcation lesions. Neth Heart J 2012 Jul 5. [Epub ahead of print]

Acute Coronary Syndrome: 42%

Baseline Characteristics	N=905
Age (years)	65±10
Male gender	76%
Diabetes Mellitus	23%
Hypertension	63%
Current smoker	20%
Previous MI	30%
Prior revascularization	40%
Indication for PCI	
Stable Angina	52%
Unstable Angina	23%
NSTEMI	13%
STEMI	6%

Complex Bifurcations: 85%


Lesion Characteristics (visual estimate)		N=929
Bifurcation Location	Left main	5%
	LAD / diagonal branch	67%
	Cx / obtuse marginal	16%
	Posterolateral / posterior descending	10%
True bifurcation	(Medina 1,1,1; 1,0,1; 0,1,1)	85%
Side branch angles	Narrow (<30°)	22%
	Large (>30°)	10%
Reference vessel diameter	Main branch	3.0 [3.0-3.5]
	Side branch	2.5 [2.5-2.7]
Lesion length	Main branch	20 [13-33]
	Side branch	7 [4-12]

Device Not Approved for Sale in the US and Japan.

www.trytonmedical.com

TRYTON[®]
Medical
BUILT FOR BIFURCATION


Tryton Clinical TVF Rates Compare Favorably with Historical Data



1. Romagnoli, et al. Real-world outcome of coronary bifurcation lesions in the drug-eluting stent era - Results from the SICI-GISE I-BIGISAm Heart J. 2010
2. Brodie et al. Outcomes and Complications With Off-Label Use of Drug-Eluting Stents: Title and Results From the STENT (Strategic Transcatheter Evaluation of New Therapies) Group JACC int. 2008
3. Hakeemet et al. Provisional vs. Complex Stenting Strategy for Coronary Bifurcation Lesions: Meta-Analysis of Randomized Trials J INVASIVE CARDIOL 2009
4. Garg et al. Implantation of the BES in patients with high SYNTAX score is associated with decreased cardiac mortality compared to a permanent polymer SES: two year follow-up results from the "all-comers" LEADERS trial - EuroIntervention 2011


TLR Rates of 2.9% at 6 Months and 4% at 12 Months

Tryton Clinical Results Challenge the Paradigm of Provisional Strategy



Low Thrombosis Rate of 0.5% at 6 and 12 Months

Tryton Thrombosis Rates Demonstrate its Safety Profile



1. Data presented at EuroPCR 2012. 2. TVF rate indicated. 3. Combined 1 and 2 stent strategy. *Data on File at Tryton Medical.

Device Not Approved for Sale in the US and Japan.

Corporate Headquarters
1000 Park Forty Plaza, Suite 325
Durham, NC 27713 USA
Phone: +1-919-226-1490
Fax: +1-919-226-1497

Tryton Medical B.V.
Centaurusweg 123
5015 TC Tilburg,
The Netherlands

To learn more visit: www.trytonmedical.com

TRYTON[®]
Medical

BUILT FOR BIFURCATION

©2012 Tryton Medical, Inc. All rights reserved. D2091G_V01